



## UDS-P LED

- Detachable handpiece with light, better visibility.
- Scaling, perio, endo.
- Automatic frequency tracking.
- Accurate feedback control and constant power output, ensure a more comfortable feeling to patients.

|                  |                              |                       |
|------------------|------------------------------|-----------------------|
| Power input      | 110V~ 50Hz/60Hz              | 220 - 240V~ 50Hz/60Hz |
| Main unit input  | 24V~ 50Hz/60Hz               | 1.3A                  |
| Output power     | 3W~20W                       |                       |
| Frequency        | 28kHz±3kHz                   |                       |
| Water pressure   | 0.1bar~5bar (0.01MPa~0.5MPa) |                       |
| Main unit weight | 0.68kg                       |                       |
| Adapter weight   | 1.2kg                        |                       |
| Dimensions       | 178mm×145mm×89mm             |                       |

### Scaling

|                       |  |  |  |
|-----------------------|--|--|--|
| G1                    |  |  |  |
| Supragingival scaling |  |  |  |
| G2                    |  |  |  |
| Supragingival scaling |  |  |  |
| G3                    |  |  |  |
| Supragingival scaling |  |  |  |
| G4                    |  |  |  |
| Supragingival scaling |  |  |  |
| G5                    |  |  |  |
| Supragingival scaling |  |  |  |
| G6                    |  |  |  |
| Supragingival scaling |  |  |  |

### Periodontics

|                     |  |  |  |
|---------------------|--|--|--|
| P1                  |  |  |  |
| Periodontal scaling |  |  |  |
| P3                  |  |  |  |
| Periodontal scaling |  |  |  |



|  |  |  |  |
|--|--|--|--|
| G7   |  |  |  |
| Remove dental crown  |  |  |  |
| G8   |  |  |  |
| Remove dental pile, with irrigation  |  |  |  |
| G9   |  |  |  |
| Supragingival scaling  |  |  |  |
| G10  |  |  |  |
| Supragingival scaling  |  |  |  |
| G11  |  |  |  |
| Orthodontic treatment  |  |  |  |
| P4   |  |  |  |
| Periodontal scaling  |  |  |  |
| P3D  |  |  |  |
| Smoothing restoration overhangs and debride in-depth root (diamond-coated) |  |  |  |